ΔVANOS*

FARRELL* VALVE
CLOSED ENTERAL
DECOMPRESSION SYSTEM
WITH ENFit® CONNECTORS

Instructions for Use

AVANOS* FARRELL* Valve Closed Enteral Decompression System with ENFit® Connectors


Rx Only: Federal Law (USA) restricts this device to sale by or on the order of a physician.

- ENFit® design employed for reducing the risk of tubing misconnections
- For use with enteral feeding devices
- FOR ENTERAL USE ONLY, do not administer IV medications through this enteral device
- The FARRELL* Bag is vented to allow the escape of gas, while also providing a container to retain excess formula in a closed system
- FARRELL* Tubing is marked in 1/4 ml increments
- Change in conjunction with enteral administration set – maximum use 24 hours


Indications for Use

The AVANOS* FARRELL* Valve Closed Enteral Decompression System with ENFit® Connectors is a closed reservoir system intended to allow the evacuation of excess gas (gastric distention/bloating) and drainage/collection of enteral feeding formula and gastrointestinal tract contents from patients utilizing an enteral device. The FARRELL* Valve System with ENFit® Connectors is for use with neonates, pediatric, and adult patients.


1. Hang FARRELL* Bag at the same height as the feeding bag/container.
2. Close WHITE clamp above FARRELL* "Y" Port.


3. a. Attach enteral administration set connector to "Y" Port on FARRELL* Tubing.
Note: For syringe feed, attach syringe to "Y" Port.
- b. Prime FARRELL* Tubing to line (3 cm) before distal end of FARRELL* Set Connector.
Note: For continuous feeding, prime pump according to pump instructions.
- c. Close BLUE clamp.

- Attach FARRELL* Set Connector to feeding tube.
(Reminder: Confirm connection to an enteral feeding device only and NOT to an IV set.)
- Important:** Position the FARRELL* “Y” port at or below the patient’s stomach.

- Open the BLUE clamp to establish flow, then open the WHITE clamp.
Note: Normal height of formula in the FARRELL* Tubing will be slightly above patient’s stomach level. Formula may continuously move up and down in the FARRELL* Tubing.

Refer to Tip #4 on the next page for additional information on establishing flow.

⚠ Warning: When FARRELL* Tubing line is open (Open WHITE clamp), the pump’s “occlude” alarm will not function, as formula will continue to flow into FARRELL* Bag.

Medication Administration

- When administering medication, use the access port on the feeding tube if possible.
- Close the BLUE clamp before opening the feeding tube access port.
- Administer medication.
- Wait 5–10 minutes before closing the feeding tube access port and re-opening the BLUE clamp.

Total FARRELL* Line Volume	9.6 ml
Priming Volume**	3.2 ml

**Volume from “Y” Port to FARRELL* Set Connector

Additional Tips and Information

1. In order for the FARRELL* Valve System to work properly, the FARRELL* “Y” port must be AT or BELOW the patient’s stomach.
2. If a patient has a distended stomach, manually decompress the patient’s stomach with a syringe prior to the initial use of the FARRELL* Valve System.
3. It may be helpful to prime the FARRELL* Tubing slightly above the in-line “Y” Port to the expected level of the fluid in the patient’s stomach.
4. Patients on very low volume feeds (for example: neonates) may require substantially longer time (60–90 minutes) to establish flow before opening the WHITE clamp above the “Y” Port.
5. When the FARRELL* Valve is working properly, formula may continuously move up and down in the FARRELL* Tubing.
6. The FARRELL* Valve System continuously decompresses the patient’s stomach utilizing a vent to the ambient environment while also providing a reservoir bag to retain excess fluid (formula, electrolytes, and/or medications) in a closed system. After relief of gastric pressure, the collected fluid can be delivered to the patient to ensure delivery of prescribed nutrients and medications.

Note: *The FARRELL* Bag will not usually inflate with gas because it is vented.*

7. If a feeding tube occlusion is suspected, close the BLUE clamp and disconnect the FARRELL* Valve System. Check for occlusions in the feeding tube. Once the occlusion is resolved reconnect the FARRELL* Valve System.
8. Discard the FARRELL* Valve System (including any residual fluid) with the enteral administration set (maximum use 24 hours).

Important Note: *If a large volume of fluid is in the FARRELL* Tubing and Bag after resolving an occlusion, adjust the head height of the FARRELL* Bag to prevent a free-flow bolus of fluid into the patient.*

⚠ Caution: While the FARRELL* Valve is in operation, flow to the patient is essentially controlled by gravity. With normal (Non-FARRELL*) operation, pump pressure (12–18 psi/ 83–124 kPa) may overcome and clear occlusions or kinks in the feeding tube. With the FARRELL* Valve in operation, these occlusions may not be overcome and formula will back up into the FARRELL* Bag. If formula begins to back up into the FARRELL* Bag, close the WHITE clamp located on the line above the FARRELL* “Y” Port and wait several minutes to determine if the pump can then overcome the occlusion. Make sure the WHITE clamp above the FARRELL* “Y” Port is CLOSED and the BLUE clamp below the FARRELL* “Y” Port is OPEN.

Interconnectability to Non-Enteral Medical Devices

ENFit® satisfies the ISO standard (80369) for small bore enteral devices. However, the design of the ENFit® Connector cannot overcome all chances of misconnection. The following connector types are potential misconnections for the ENFit® Connector (feeding set and “Y” Connectors) of this enteral feeding tube:

- Suction ports on Endotracheal Suction Systems
- Respiratory circuit filtration connectors
- Oxygen inlet connectors for Resuscitation Devices
- Baxter® IV Solution Bag ports (such as NaCl, Ringers Solution, etc.)
- Sample ports on drainage bags
- Peritoneal Dialysis connectors
- Cones & sockets of ISO 5356-1:2004 & ISO 5356-2:2004
- Temperature sensor connectors & mating ports of ISO 8185:2007
- Oxygen nipples as defined in EN 13544-2:2002
- Nebulizer system ports
- IV Spike Syringe Ports
- Bladder Drainage Catheter Outlet Ports

Note: Infrequently occurring factors that may result from the failure or misuse of a feeding decompression accessory device include: negligible tissue irritation or patient discomfort, contamination, tissue necrosis, fluid aspiration, inhibited lung function, medication/nutrition misdose or delay and related complications, airway obstruction, or allergic reaction.

 Single Patient Use	 Do not use if package damaged	 Product is not made with DEHP as a Plasticizer
 Product is not made with BPA		 Not made with natural rubber latex
RX Only	 Caution	 Consult instructions for use

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www.avanos.com

Distributed in the USA by Avanos Medical Sales, LLC,
5405 Windward Parkway, Alpharetta, GA 30004 USA.
In USA, 1-844-4AVANOS (1-844-428-2667). www.avanos.com

 Avanos Medical, Inc., 5405 Windward Parkway, Alpharetta, GA 30004 USA.

 Avanos Medical Belgium BVBA, Leonardo Da Vincilaan 1, 1930 Zaventem, Belgium.

Sponsored in Australia by Avanos Medical Australia Pty Ltd,
475 Victoria Avenue, Chatswood, NSW 2067, Australia.

 2797

*Registered Trademark or Trademark of Avanos Medical, Inc., or its affiliates.

♦Baxter is a Registered Trademark or Trademark of Baxter International, Inc.

© 2018 AVNS. All rights reserved. 2020-04-27 15-M1-325-02 / 70215483